

RNA Sequencing

Elisha Roberson, Ph.D.


Depts. of Medicine & Genetics

eroberson@wustl.edu

2016-December-06


Comparing RNA-Seq and microarrays

RNA-Seq vs. microarrays




- Hybridization requires known targets

RNA-Seq vs. microarrays




- Expression detected by fluorescence

RNA-Seq vs. microarrays




- Low-level expression can be difficult to detect compared to background

RNA-Seq vs. microarrays



- High-level expression can saturate probes

RNA-Seq vs. microarrays



- Variation in the subject's RNA sequence can affect binding kinetics

RNA-Seq vs. microarrays

RNA-Seq differences

- Counting – no probes to saturate (though some sequences can predominate the library) and can always sequence more to get low-expressors
- Doesn't *require* a known gene sequence
- Not affected by variation, as long as it doesn't affect transcript stability

Isoform A



Isoform B




Sample prep is ***critical***

Collection & storage



- Blood is low-risk, easy to access, and frequently used
- ***But*** strongly enriched for red cells. What transcripts will predominate?

Collection & storage



- Freeze PBMCs directly
- Lyse directly in Trizol


Collection & storage




Tissues

- Biopsies, whole organs, etc.
- Unlike blood, the complex structure and embedded connective fibers make tissues more difficult to process and preserve.

Collection & storage




- Some solutions, like RNAlater, can preserve tissue RNA for short-term storage at RT or long-term storage at -20°C



- PBMCs shipped on (too little) dry ice
- Degraded 3'-5'
- PolyA kit fails


Library prep / cDNA methods

Polyadenylation (polyA) preps



- polyT primed first-strand synthesis works
- **But** can lead to 3' bias and only captures polyadenylated transcripts.


Polyadenylation (polyA) preps




Wash gently




Pipette vigorously



Random hexamers




cDNA




Ribosomal depletion preps


- Majority of cellular RNA is non-coding, particularly ribosomal RNA
- RNA polymerase I (28S, 18S, 5.8S rRNA) and Pol III (5S rRNA)
- Not polyadenylated



Ribosomal depletion preps




Ribosomal depletion preps




Stranded preps

- Standard cDNA → library prep retains no information about transcript strand.
- Some loci have antisense transcripts



Stranded preps



- Normal first strand synthesis. 2nd strand incorporates uracil
- Uracil-DNA glycosylase excises U-base from DNA
- Endonuclease VIII breaks backbone at those sites

Sequencing choices

Short reads

- Illumina
 - Single-end vs. paired-end
 - Paired-end superior. Estimates insert size empirically
 - Read length
 - Greater cycle number preferred
 - 2x75 good compromise
 - Depth / coverage
 - Very different from DNA seq
 - Variable gene length and expression level
 - Several tools to estimate

Long reads – Pac Bio SMRT

- Full-length isoform sequencing (\$\$\$\$\$)

<i>Pacbio Library Construction and Sequencing</i>	<u>Cost Per Sample</u>	
Sequencing SMRT Cell	\$257	Min \$1,100 / sample for whole transcriptome.
Standard Library Prep	\$560	
Low_input Library Prep	\$603	
Iso-Seq Whole Transcriptome Lib_Prep	\$875	
Iso-Seq Targeted Lib_Prep	\$664	

Long reads – Oxford nanopore

- MinION / PromethION (also \$\$\$\$\$\$)



FLO-MIN106

SpotON Flow Cell Mk I (R9.4)

\$900.00

Single

12-pack

24-pack

48-pack

— 1 +

Buy now



SQK-LSK208

Ligation Sequencing kit 2D (R9.4)

\$599.00

— 1 +

Buy now

[Product overview](#)

The Ligation Sequencing Kit 2D is designed to prepare genomic, amplicon and cDNA, with or without barcoding, for sequencing on the Oxford


Alignment and counting

RNA-Seq does not look like the genome

Genome locus




Isoform A




Isoform B




RNA-Seq does not look like the genome



RNA-Seq does not look like the genome



RNA-Seq does not look like the genome



Genomic aligners expect the library to reflect **genome** architecture. Intron splicing looks like large deletions, and can confuse aligner.

One alternative is to align to transcript FASTA rather than whole-genome.

Splice-aware genome aligners

- Tophat2 (bowtie derived)
- Spliced Transcript Alignment to Reference (STAR)


Transcript alignment

- Kallisto
- Sailfish / Salmon

De novo assembly

- ABySS / TransABySS
- Trinity
- SOAPdenovo-Trans








Ribosomal depletion alignment



Counting reads

Counting tools


HTSeq count
subread featureCount
RSEM

	union	intersection_strict	intersection_nonempty
	gene_A	gene_A	gene_A
	gene_A	no_feature	gene_A
	gene_A	no_feature	gene_A
	gene_A	gene_A	gene_A
	gene_A	gene_A	gene_A
	ambiguous	gene_A	gene_A
	ambiguous	ambiguous	ambiguous

*HTSeq count


RNA-Seq complications - duplicates

Genome locus




Highly covered genomic locus for DNA-Seq


RNA-Seq complications - duplicates



RNA-Seq complications - duplicates



RNA-Seq complications - duplicates




RNA-Seq complications - duplicates



- For DNA-Seq the target is equimolar, but RNA-Seq is more complicated
- Important considerations:
 - Sequencing depth
 - Gene relative expression level
 - Gene size

RNA-Seq complications - duplicates



- **But** this is a solvable problem
- Adding a short barcode to each fragment during PCR, called a unique molecular identifier (UMI) we know whether reads are truly unique, even with identical 5' mapping

Estimating gene abundance & differential expression

Abundance metrics


- Fragments per Kb of exon per million mapped reads (FPKM)
- Transcripts per million (TPM)
- TPM preferable
 - Different total reads between experiments skew FPKM
 - TPM consistent, i.e. 1 TPM sample A and sample B really means similar abundance

Modeling counts – edgeR, DESeq2

- Counts are not normally distributed
- What models counts?
 - Poisson distribution
 - But assumes variance & mean equal
 - Negative binomial
 - Mean \neq variance
 - edgeR, DESeq2 R packages

Modeling counts – cuffdiff2


4) Combine uncertainty and overdispersion into a single model of fragment count variability (beta negative binomial)



5) Test for significance of changes between conditions in transcript-level counts


- edgeR & DESeq2 model gene-level differential expression
- cuffdiff2 tests for significant isoform-level DE

Modeling counts - VOOM




- VOOM uses \log_2 of counts per million normalization factor
- Differential expression using the empirical Bayes limma pipeline

How similar are gene DE algorithms?




Transcript k-mer modeling - sailfish




- Hashing uses a large amount of memory
- But lookups are blazingly fast
- Calculates TPMs

Transcript pseudoalignment - kallisto




- Builds a de Bruijn graph of transcript sequences
- Pseudoalignment – compatible transcripts, not where in transcript
- Very fast and efficient
- Allows for bootstrapping

Transcript DE from kallisto - Sleuth




Single cell RNASeq – Fluidigm C1





- Microfluidics capture single cells
- Lyse cells and generate cDNA
- Requires live cells
- 96-800 cells / chip

Single cell RNASeq – DropSeq




- Flows beads in a droplet.
- Cells, usu singles, merge into a droplet for library prep.
- 10s of thousands of cells.

RNA-Seq to diagnose Mendelian disease



c)



Summary

- RNA-Seq vs. microarrays
 - Microarray requires knowing target sequence
 - Poor dynamic range (hard to detect low-level expressors, saturate at high-levels)
- Sample collection & storage
 - Blood should be spun down to PBMCs
 - Can be directly lysed
 - Tissue should have RNAlater applied soon as possible, followed by disruption.

Summary

- cDNA / library prep
 - Polyadenylation library methods selectively capture polyA transcripts
 - Ribosomal depletion methods degrade ribosomal RNA, but leave non-polyadenylated
 - Strandedness
 - There are an appreciable number of genes with antisense transcripts.
 - Also useful for identifying genes in species without a reference genome

Summary

- Sequencing technologies
 - Long reads (expensive), but sequence full isoform
 - Short read. Reasonable price.
- Aligner
 - Must use an aligner that is aware of introns
 - May align to either genes (STAR, etc) or transcriptome (Tophat 2 and kallisto)
- Counting
 - Subread featureCounts, HTSeq count, RSEM

Summary

- Differential expression
 - Gene, negative binomial: edgeR, DESeq2
 - Gene, log2 counts per million: VOOM
 - Transcript, TPM, kallisto.
 - Transcript, TPM, Sailfish
- Single-cell
 - Isolate a single cell and make a libraries (low-output)
 - Spike-ins help
 - UMIs help also